

Mesocosm User Manual

Dr. Nyssa Silbiger and Danielle Barnas

2019-05-03

Contents

1	Summary	5
2	Contacts	7
3	System Details	11
4	Inventory	15
5	Start-up Guide	17
6	Tidal Manipulation	19
7	Controlling pH	21
8	Apex Programming Guide	23
9	Apex Fusion	27
10	Breaker Box Connections	29

Chapter 1

Summary

This manual describes the design, operation, and maintenance of the mesocosm aquaria, located in the loading bay between Citrus Hall and Eucalyptus Hall at California State University, Northridge, funded and operated by Dr. Nyssa Silbiger.

Chapter 2

Contacts

Name	Involvement	Contact	
		Info	Notes
Dr. Nyssa Silbiger	System De- sign Asst. Pro- fessor, CSUN	nyssa. silbiger@ csun. edu 818- 677- 4427	
Danielle Barnas	System Instal- lation and Main- te- nance Sil- biger Lab Tech, CSUN	danielle. barnas@ csun. edu	
Louis Dang	Systems Engineer	louis@ aqualogicinc. com	www.aqualogicinc.com
Bill Krohmer	Administra- tion Operations	william. krohmer@ csun. edu	

Name	Involvement	Contact	
		Info	Notes
Science Shop	CSUN College of Science and Math Machine Shop	818-677-3055	Location: EH 2014 Available M-Th 0600-1630 www.csun.edu/Science-Shop
Perry	Supervising	perry.	
Martin	Plumber	martin@PPM	
		csun.edu	
		818-677-2222	
Will	Network	will.	
Moran	Engineer	moran@csun.	
	CSUN IT	edu	
		818-677-6273	
Willy	Lead	willy.	
Martinez	Electrician	martinez@csun.	
	PPM	edu	
	Electric Shop	818-677-6273	
Neptune Systems	Apex Support Team	support@www.neptunesystems.com	
		neptunesystems.com	

Contact			
Name	Involvement	Info	Notes
Dickson	Seawater	co2crms@	Marine
Lab	CO2	ucsd.	Physi-
	CRMs	edu	cal
		858-	Lab
		534-	Scripps
		2582	Insti-
			tution
			of
			Oceanog-
			raphy
			Uni-
			ver-
			sity of
			CA,
			San
			Diego
			9500
			Gilman
			Drive
			La
			Jolla,
			CA
			92093-
			0244
			USA

System Details

Aquaria System Component List

- Experimental Tanks 21.25x12.5x13.5H - Per tank:
- 1 Submersible powerhead pump [*HydorNanoKoralia*240powerhead](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed8765e4152fd9dd11dc)
 - 1 200 W Heater [*Hidoraquariumheater*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed8765e4152fd9dd11dc)
 - 1 Light [*HaloBasicM-110*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed8765e4152fd9dd11dc)
 - 1 Temperature probe *Apex*
 - 1 pH probe *Apex*
 - 1 Solenoid valve for pH *Apex*
 - 3 Flow sensors *Apex*, FS – 251/4fitting, flowratesfrom3 – 12gph
 - 1 Main Supply line: “N”
 - 1 Solenoid Supply line: “S”
 - 1 Drain line: “D”
 - 1 Gate valve solenoidforwaterinflow
 - 1 VDM [*ApexVariableDimmingModule*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed8765e4152fd9dd11dc)
 - 1 FMM [*ApexFluidMeteringModule*](<https://www.neptunesystems.com/getstarted/fmk/>)
 - 1 PM1 [*ApexProbeModule1*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed8765e4152fd9dd11dc)
 - 1 Base Unit [*Apex*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed8765e4152fd9dd11dc)

- 1 EB832 [*Apex8-outletEnergyBar*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed870)
- 1 CO2 regulator valve [*TunzePHControllerSet, pressure reducing valve7077/3*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed870)
- 1 Industrial Grade Carbon Dioxide, 50 pound cylinder

Filtration and Recirculation System Component List

- Sump 66.25x31.5x21H
- Chiller [*AquaLogicMulti-TempandTitanSeries*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed870)
- Heat Pump [*AquaLogicMulti-TempandTitanSeries*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed870)
- PM1 *ApexProbeModule1*
- PhosBan chemical filter [*PhosBanReactor550*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed870)
- Water pump [*PerformanceProCascadepump*](https://github.com/SilbigerLab/Mesocosm_UserManual/tree/394a3f7d9fed870)
- UV Sterilizer *CometSeries95WattLamp*
- CO2 Scrubber 4outflowtubingports
- Airstones 4units
- Carbon filter cells 3units, CF28AC, 28in, ActC
- Mesh filter 8units, MatalaFilterMedia, interchanged4Bluehighdensityand4Blacklowdensitysheets

System Operational Sequence

- This is a closed loop system where water from each individual tank will recirculate back to a main holding reservoir *sump*.
- Normal High Tide operating water level is approximately 12.5"H for a total water volume of 14.37 gal per tank 287.4galtotalforthe20 – tank – system.
- Normal Low Tide operating water level is approximately 4"H for a total water volume of 4.60 gal per tank 92.0galtotalforthe20 – tank – system.
- Excess water volume to sump at low tide is 9.77 gal per tank 195.4galtotalforthe20 – tank – system.
- Normal sump operating water level is 7" water in the filter cell compartment, which has an approximate volume of 82.32 gal. Sump freeboard volume is 107.39 gal.
- Aquaria drain line is in line with a 30 gal sump pump, which will draw water from the aquaria drain line and pump water to the sump.
- Sump is in line with a secondary holding tank for sump overflow at Low Tide.
- 195.4 gal returning to sump in a Low Tide scenario
- 30 gal in sump pump
- 107.39 gal in sump *freeboardvolume*
- 58.01 gal necessarily pumped to secondary holding tank
- Flow rate for each tank is 2-6 GPH.

- Chiller is plumbed inline with the filtration skid which includes mechanical/biological filtration as well as UV sterilization *chemical filtration*.
- One main pump recirculates the water flow throughout the experimental tanks and the main holding reservoir.
- Each tank has an immersion heater that allows tank temperatures to be set 15 degF 8.3degC above the main holding tank reservoir.
- **Note: The tank needs to be static in order to heat up to a desired temp. Once the temperature has been reached then it can be set to flow through mode.**
- A small submersible powerhead in each tank provides water circulation throughout the tank.
- Each tank has 2 supply lines, each with 1 flow sensor, and 1 gate valve in line with 1 supply line for tidal effect. Each tank also has 2 drain lines with 1 flow sensor. Incoming and outgoing flow rates have to be manually adjusted for the tidal effect.
- outgoing tide: incoming flow rate to be lower than outgoing flow rate.
- incoming tide: incoming flow rate to be greater than outgoing flow rate.
- Note: The tank will not be completely empty during low tide events to prevent the recirculating powerhead from running dry.
- Each tank has 1 CO2 supply line with an airstone connected to 1 solenoid valve to lower pH in tanks.
- Air compressor connected to a CO2 scrubber will bubble air into the sump to bring pH to ambient or near-ambient conditions in the holding reservoir.
- Each tank has individual LED lighting.
- Tanks controlled via Neptune Systems Apex Controllers. Each Apex controls 4 tanks.
- Controllable parameters are pH, temperature, tidal effect, and lighting.

Apex Connection Series

- Each EnergyBar connects to the Base Unit with an AquaBus cable via the AquaBus Ports for power. 2 EB832 units connect to 1 Base Unit.
- Each CO2 Solenoid valve connects to the EnergyBar via the DC24 Accessory Port on the side of the EB832. 2 Solenoid valves connect in 1 EB832.
- 1 PM1 connects to 1 EnergyBar with an AquaBus cable via the AquaBus Ports, and all PM1 modules connect in series with each other for power.
- VDM connects to the last PM1 in series with an AquaBus cable via the AquaBus Ports for power.
- Temperature probes connect to the PM1 Temp Port or the Base Unit Temp Port. 1 Temperature probe in each PM1, and 1 Temperature probe in the Base Unit.
- pH probes connect to the PM1 pH/ORP Port or the Base Unit pH/ORP Port. Push the BNC female connector of the probe on to the male connector and turn 1/4 turn clockwise to lock the connector in

place. 1 pH probe in each PM1, and 1 pH probe in the Base Unit.

- Halo light connects to the VDM or Base Unit via the V1/V2 or V3/V4 Port. 2 Light connections in the VDM and 2 Light connections in the Base Unit.
- 1 FMM connects to 1 EnergyBar *whichever EB832 is not powering the PM1 modules* with an AquaBus cable via the AquaBus Ports, and all FMM connect in series with each other for power.
- 3 Flow sensors connect to each FMM via 3 of the numbered ports.

EB832 Outlet Connections

Note: Each horizontal row on an EB832 corresponds to one tank, yielding 4 outlets per aquarium. In order:

1. 200W Heater
2. Hydor Powerhead
3. Water supply line “S” Solenoid
4. Halo Light

Chapter 4

Inventory

Experimental Mesocosm

Item	Quantity
Experimental Aquarium	20
Aquarium Lid	20
200W Heater	20
Suction Heater Slip	40
Halo Light	20
Halo Light Power Cords	20
Hydor Powerhead	20
Temperature Probe	20
Suction Temperature Probe Slip	20
pH Probe	20
Suction pH Probe Slip	20
pH Calibration Pack 4, 7, 10	10
Solenoid Valve	20
Solenoid Power Supply	20
Flow Sensor	60
Gate Valve	20
Inflow Seawater Tubing	40
Inflow CO2 Tubing	20
Airstone	20
Inflow Tubing Stand	40
VDM	5
FMM	20
PM1	16
Apex Base Unit	5
Apex EB832	10
Tunze CO2 Regulator Valve	1
180gal Sump	1
Chiller	1
Heat Pump	1

Filtration

Item	Quantity
PhosBan Chemical Filter with Media	1
Cascade Water Pump	1
UV Sterilizer	1
UV Light	1
CO2 Scrubber	1
Airflow Tubing	4
Airstone	4
Carbon Filter Cell	3
Matala Filter, High Density <i>Blue</i>	4
Matala Filter, Low Density <i>Black</i>	4

Spare Items

Item	Quantity
Solenoid Valve	4
Airstone	29
Halo Light Cable	16

Chapter 5

Start-up Guide

1. Operating water level in the filtration sump should be 7" in the filter cell compartment.
2. Water from the container will feed back to the filtration skid, and from there it is pumped through the UV sterilizer and chiller barrel, then back to the container.
3. There are two valves located after the chiller
 - a. One controls the flow directly back to the tanks
 - b. The other valve is the bypass valve which diverts the flow back to the sump This is used to regulate the line pressure going back to the container. The more closed, the higher the pressure in the line, and the more open, the lower the pressure. The chiller has a safety flow switch that requires a minimum flow rate for the chiller to operate, so the bypass valve is used in this case to regulate the chiller flow as well as the container flow.
4. Before filling the tanks make sure the drain valve located under the tank is closed.
5. Fill each rack one at a time and make sure rack and filtration skid flows are balanced before moving on to the next rack.
6. Make sure the complete system reaches equilibrium in standard recirculation mode before setting up the tidal cycle.

Chapter 6

Tidal Manipulation

Controlling the tidal cycle of each experimental tank with the Apex. This is achieved by manipulating the incoming and outgoing flow rates of each individual tank with the needles described in the

SystemDetails

, and setting the ON/OFF time cycle of the supply line with the solenoid. The basic procedure is outlined below.

1. Set the flow rate of the supply line N

#

FLW, the one without the solenoid, for example 12 Liters/Hr, by slowly turning the black knob.

- a. Note that the Apex controller has some lag time in registering the flow rate after the valve has been adjusted, the delay can be up to 30 seconds or more. Once the rate is set you should check periodically to make sure the rate has not changed using a graduated cylinder and a timer.

2. Adjust the outgoing flow rate of the drain line D

#

FLW higher than the N

#

FLW, for example 19.04 Liters/Hr.

- a. With the above condition, the outgoing flow rate is higher than the incoming, so this will create the low tide effect.

3. To set the high tide effect, manually turn on and adjust the flow rate of the supply line S

#

FLW , for example 14.08 Liters/Hr.

4. Once the S

#

FLW is set, change setting of SOL-TNK-# outlets3and7oneachEB832 to AUTO on the Fusion page or using the Display Module. For a constant ON/OFF cycle over a 12.5 hour period, the Advanced program should look something like the program below.

Fallback ON

Osc 000:00/375:00/375:00 then ON

1. In the event the EnergyBar loses connection with the Apex Base, Fallback: ON will keep the solenoid open, allowing water to continuously flow from S

#

FLW.

2. The oscillate command as written will open flow from S

#

FLW for 6.25 hours, initiating the High Tide scenario, then close for 6.25 hours, initiating the Low Tide scenario. This will provide the effect of two high and two low tides of a semidiurnal tidal cycle over a 25 hour period.

3. Using the flow rates stated above, each tidal shift will last 5.25 hours and maintain the tide for 1 hour.

For more advanced programming features, see the Comprehensive Manual. Start on Page 65 for Seasonal Features and Moon cycles.

Chapter 7

Controlling pH

The pH is controlled with the addition of CO₂ gas to the system. The gas is delivered to the tank by air stone and is controlled through the Apex Controls with a solenoid valve connected to the EB832.

1. Once the CO₂ regulator is connected to a tank, open the main tank valve.
2. Use the pressure adjusting screw to adjust the pressure *inbar* on the pressure gauge. Turning **clockwise to open**, thus increasing pressure, while turning **counterclockwise to close**, thus reducing pressure.
3. The pressure should be set to 0.5 up to 1 bar on the gauge *7.5psi*
4. Open the fine adjustment valve to allow gas to the tank solenoid. If the pressure on the gauge is too high this may prevent the CO₂ solenoid from completely closing, which will inject excess CO₂ into the system.
5. Programming the solenoid: pH-TNK-#

Control type: pH Control

Probe name: pH

Fallback: OFF High Value: 8.2

Low Value: 7.9

On when: High

Refer to Comprehensive Manual for set point programming.

Chapter 8

Apex Programming Guide

Recommendations for programming the Apex aquarium controllers designated for the Silbiger Lab Mesocosm, located in the loading bay between Citrus Hall and Eucalyptus Hall at California State University, Northridge. These recommendations are for maintaining tanks at ambient conditions. Changes should be made according to your study aims.

The following are using the numbered system of Apex_39106, controlling tanks 1-4. All methods are transferrable across all 5 Apex controllers to yield the same outcome in all 20 tanks.

Contents

Probes

Outlets and Ports

Outlet Setup in ApexFusion

Profiles

Probes

- Tmp-1 *Base*
- pH-1 *Base*
- Tmp-2 *PM1_2*
- pH-2 *PM1_2*
- Tmp-3 *PM1_3*
- pH-3 *PM1_3*
- Tmp-4 *PM1_4*
- pH-4 *PM1_4*

Outlets and Ports

- Base Unit
- WHITE-TNK-1
- BLUE-TNK-1
- WHITE-TNK-2
- BLUE-TNK-2
- SndAlm_I6
- SndWrn_I7
- EmailAlm_I5
- Email2Alm_I9
- EB832_1
- HEATER-1
- PWRHD-1

- SOL-TNK-1
- LIGHT-TNK-1
- PH-TNK-1
- HEATER-3
- PWRHD-3
- SOL-TNK-3
- LIGHT-TNK-3
- PH-TNK-3
- EB832_2
- HEATER-2
- PWRHD-2
- SOL-TNK-2
- LIGHT-TNK-2
- PH-TNK-2
- HEATER-4
- PWRHD-4
- SOL-TNK-4
- LIGHT-TNK-4
- PH-TNK-4
- VDM
- WHITE-TNK-3
- BLUE-TNK-3
- WHITE-TNK-4
- BLUE-TNK-4
- BluLED_11_5
- WhtLED_11_6
- FMM_1
- S1-FLW
- N1-FLW
- D1-FLW
- FMM_2
- S2-FLW
- N2-FLW
- D2-FLW
- FMM_3
- S3-FLW
- N3-FLW
- D3-FLW
- FMM_4
- S4-FLW
- N4-FLW
- D4-FLW

Outlet Setup in ApexFusion

- LIGHT-TNK-#
- Fallback OFF
- Set OFF
- If Sun 0/0 Then ON
- If Moon 0/0 Then ON
- WHITE-TNK-#
- Fallback OFF
- Set OFF
- If Sun 0/0 Then RampUp

- BLUE-TNK-#
- Fallback OFF
- Set OFF
- If Moon 0/0 Then RampUp
- WhtLED_#
- Fallback OFF
- Set OFF
- If Sun 0/0 Then RampUp
- BluLED_#
- Fallback OFF
- Set OFF
- If Moon 0/0 Then RampUp

Profiles

- RampUp:
- Ramp Time: 30 min
- Start Intensity: 0
- End Intensity: 100

Chapter 9

Apex Fusion

To access the Silbiger Lab Fusion account, click [here](#), click “Get Control” and enter the login information:

Username: SilbigerLab

Password: silbigerlab

Contents

Dashboard

Outlet Setup

Downloading Data Logs

Dashboard

- On your Dashboard you will see every outlet available to you and the current state of that outlet *ON* or *OFF*, as well as the current readings for any probes enabled on the Apex.
- On the left are graphical logs of Salinity, pH and Temperature in the first of all four tanks controlled by the Apex
- The top menu bar options *left to right*:
- Apex List: returns you to the menu of available Apexes linked to your Fusion account
- Alarm Log: shows any daily alarms triggered by the system, as set by your parameters
- Input Log:
-

Outlet Setup

1. Click Expand along the top bar *depicted as three gears*.
2. Immediately next to the gear symbol is the Outputs icon *depicted as a three pronged outlet*. Click this icon to view your array of outlets and other outputs from the system, as well as manipulate these items. These are arranged by the name you give them, the device they’re connected to, the type of output, and whether or not you have chosen to log this item’s history.
3. Click any “outlet” type to

Downloading Data Logs

Format of what to enter into your internet browser:

<http://:/cgi-bin/outlog.xml?sdate=yymmddhhmm&days=n>

Examples of what to enter into your internet browser:

<http://130.166.116.174/cgi-bin/outlog.xml?sdate=190426>

<http://130.166.116.174/cgi-bin/datalog.xml?sdate=190426&days=7>

- the value after sdate= is the start date for when you want logged information, and days=n yields data n days after that start date.

Chapter 10

Breaker Box Connections

Following switches from top-down, then left-right.

- 1,3: Main Disconnect
- 5,7: Air Conditioner
- 9,11: Condenser (Chiller, Filtration)
- 2: General Power (Lights and Outlet Box #1)
- 4: Outlet Box #2 (Tanks 17-20)
- 6: Outlet Box #3 (Tanks 13-16)
- 8: Outlet Box #4 (Tanks 9-12)
- 10: Outlet Box #5 (Tanks 5-8)
- 12: Outlet Box #6 (Tanks 1-4)

Powering the container: 1. Turn on (flip left to right) the Main Disconnect switches. Once powered, the remaining switches will supply power to their individual breakers. 1. Turn on (flip right to left) the General Power switch. Once powered, the light switch to the right of the entrance will turn on/off the overhead light, the O2 sensor above the light switch will be activated, and all outlets in Outlet Box #1 will be active. 1. Turn on (flip left to right) the Air Conditioner switches. Once powered, the A/C unit can be controlled via remote control or the front display panel on the unit. 1. Switches 4, 6, 8, 10, and 12 all correspond to outlet boxes lining the upper perimeter of the container. Each box is used to power up to two (2) Apex units and enough modules and other devices for up to four (4) tanks. Each outlet covering has a number corresponding to the labeling for these switches. To turn any or all on, flip the switch(es) right to left. 1. Leave the Condenser switches (9 and 11) in the “off” position unless the chiller and filtration system become connected